

Catalog #: 5127

Aliases

PDLIM5; PDZ And LIM Domain 5; LIM; Enh; Enigma-Like PDZ And LIM Domains Protein; PDZ And LIM Domain Protein 5; Enigma Homolog; Enigma-Like LIM Domain Protein; ENH1; ENH; L9

Background

Gene Name: PDLIM5

NCBI Gene Entry: [10611](#)

UniProt Entry: [Q96HC4](#)

Application Information

Molecular Weight: Predicted, 64 kDa; observed, 64 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15850

Species Reactivity: Human, rat

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human PDLIM5

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:400-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-PDLIM5 antibody (Cat#5127). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-PDLIM5 antibody (Cat#5127, 1:2,000) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using NaQ™ ECL Substrate Kit (Cat#716).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM