

Catalog #: 51278

Aliases

Inter-alpha-trypsin inhibitor heavy chain H3; ITI heavy chain H3; ITI-HC3; Inter-alpha-inhibitor heavy chain 3; Serum-derived hyaluronan-associated protein; SHAP

Background

Gene Name: ITIH3

NCBI Gene Entry: [3699](#)

UniProt Entry: [Q06033](#)

Application Information

Molecular Weight: Predicted, 99 kDa; observed, 95 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, rabbit

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human ITIH3

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of ITIH3 expression in HepG2 (A), A549 (B), rat liver (C) whole cell lysates.