

Catalog #: 51308

Aliases

RNS4; Ribonuclease 4; RNase 4

Background

Gene Name: RNASE4

NCBI Gene Entry: [6038](#)

UniProt Entry: [P34096](#)

Application Information

Molecular Weight: Predicted, 16 kDa; observed, 17 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human RNASE4

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of RNASE4 expression in PC12 (A), AML12 (B) whole cell lysates.