

Catalog #: 51349

Aliases

MCH2; Caspase-6; CASP-6; Apoptotic protease Mch-2

Background

Gene Name: CASP6

NCBI Gene Entry: [839](#)

UniProt Entry: [P55212](#)

Application Information

Molecular Weight: Predicted, 33 kDa; observed, 33,22 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Caspase 6

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of Caspase 6 expression in HEK293T (A), Hela (B), mouse brain (C), mouse kidney (D), mouse liver (E), rat kidney (F), rat liver (G) whole cell lysates.