

Catalog #: 51366

Aliases

Muscarinic acetylcholine receptor M1

Background

Gene Name: CHRM1

NCBI Gene Entry: [1128](#)

UniProt Entry: [P11229](#)

Application Information

Molecular Weight: Predicted, 51 kDa; observed, 55 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human CHRM1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of CHRM1 expression in HEK293T (A), mouse brain (B), rat brain (C), LOVO (D), MCF7 (E) whole cell lysates.