

Catalog #: 51396

Aliases

HST2; HSTF2; Fibroblast growth factor 6; FGF-6; Heparin secretory-transforming protein 2;
HST-2; HSTF-2; Heparin-binding growth factor 6; HBGF-6

Background

Gene Name: FGF6

NCBI Gene Entry: [2251](#)

UniProt Entry: [P10767](#)

Application Information

Molecular Weight: Predicted, 22 kDa; observed, 23 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, dog, pig, rabbit, sheep

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human FGF6

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of FGF6 expression in MCF7 (A), HepG2 (B), Raw264.7 (C), H9C2 (D) whole cell lysates.