

Catalog #: 51439

Aliases

FBX20; FBXO20; KIAA0858; LIM domain only protein 7; LMO-7; F-box only protein 20; LOMP

Background

Gene Name: LMO7

NCBI Gene Entry: [4008](#)

UniProt Entry: [Q8WWI1](#)

Application Information

Molecular Weight: Predicted, 192 kDa; observed, 192 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human LMO7

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of LMO7 expression in U87MG (A) whole cell lysates.