

Catalog #: 51479

Aliases

PKCB; PRKCB1; Protein kinase C beta type; PKC-B; PKC-beta

Background

Gene Name: PRKCB

NCBI Gene Entry: [5579](#)

UniProt Entry: [P05771](#)

Application Information

Molecular Weight: Predicted, 76 kDa; observed, 77,80 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, monkey, pig, rabbit

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human PKC beta

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of PKC beta expression in HEK293T (A), HeLa (B), H1688 (C), mouse lung (D), mouse spleen (E), rat lung (F) whole cell lysates.