

Catalog #: 5159

Aliases

HERC5; HECT And RLD Domain Containing E3 Ubiquitin Protein Ligase 5; CEB1; HECT Domain And RCC1-Like Domain-Containing Protein 5; E3 ISG15--Protein Ligase HERC5; Cyclin-E-Binding Protein 1; Hect Domain And RLD 5; CEBP1; Probable E3 Ubiquitin-Protein Ligase HERC5; EC 2.3.2.-; EC 6.3.2

Background

Gene Name: HERC5

NCBI Gene Entry: [51191](#)

UniProt Entry: [Q9UII4](#)

Application Information

Molecular Weight: Predicted, 117 kDa; observed, 117 kDa

Clonality: Mouse oligoclonal antibody

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human HERC5

Isotype

Mouse IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:200-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis using anti-HERC5 antibody (Cat#5159). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-HERC5 antibody (Cat#5159, 1:1,000) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM