

Catalog #: 51615

Aliases

TNT; Troponin T, slow skeletal muscle; TnTs; Slow skeletal muscle troponin T; sTnT

Background

Gene Name: TNNT1

NCBI Gene Entry: [7138](#)

UniProt Entry: [P13805](#)

Application Information

Molecular Weight: Predicted, 32 kDa; observed, 33 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, pig

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human sTnT

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of sTnT expression in H1299 (A), A549 (B), NS-1 (C), rat skeletal muscle (D) whole cell lysates.