

Catalog #: 51617

Aliases

Ubiquitin carboxyl-terminal hydrolase 44; Deubiquitinating enzyme 44; Ubiquitin thioesterase 44; Ubiquitin-specific-processing protease 44

Background

Gene Name: USP44

NCBI Gene Entry: [84101](#)

UniProt Entry: [Q9H0E7](#)

Application Information

Molecular Weight: Predicted, 81 kDa; observed, 81 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human USP44

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of USP44 expression in K562 (A), Jurkat (B) whole cell lysates.