

Catalog #: 51702

Aliases

KIAA0778; Sodium/potassium-transporting ATPase subunit alpha-2; Na(+)/K(+) ATPase alpha-2 subunit; Sodium pump subunit alpha-2

Background

Gene Name: ATP1A2
NCBI Gene Entry: [477](#)
UniProt Entry: [P50993](#)

Application Information

Molecular Weight: Predicted, 112 kDa; observed, 100 kDa
Clonality: Rabbit polyclonal antibody
Species Reactivity: Human, mouse, rat, bovine, monkey, pig
Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human ATP1A2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of ATP1A2 expression in HeLa (A), H1792 (B), Panc1 (C), AML12 (D) whole cell lysates.