

Catalog #: 51710

Aliases

JTK2; TKF; Fibroblast growth factor receptor 4; FGFR-4; CD334

Background

Gene Name: FGFR4

NCBI Gene Entry: [2264](#)

UniProt Entry: [P22455](#)

Application Information

Molecular Weight: Predicted, 87 kDa; observed, 95 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human FGFR4

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of FGFR4 expression in HepG2 (A), HEK293T (B), AML12 (C) whole cell lysates.