

Catalog #: 51712

Aliases

MLZE; Gasdermin-C; Melanoma-derived leucine zipper-containing extranuclear factor

Background

Gene Name: GSDMC

NCBI Gene Entry: [56169](#)

UniProt Entry: [Q9BYG8](#)

Application Information

Molecular Weight: Predicted, 57 kDa; observed, 58 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human MLZE

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of MLZE expression in mouse spleen (A), rat kidney (B), rat spleen (C) whole cell lysates.