

Catalog #: 51759

Aliases

GPR10; GR3; Prolactin-releasing peptide receptor; PrRP receptor; PrRPR; G-protein coupled receptor 10; hGR3

Background

Gene Name: PRLHR

NCBI Gene Entry: [2834](#)

UniProt Entry: [P49683](#)

Application Information

Molecular Weight: Predicted, 41 kDa; observed, 55,41 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human GPR10

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of GPR10 expression in HEK293T (A), Hela (B), A2780 (C), mouse brain (D), mouse kidney (E), rat kidney (F) whole cell lysates.