

Catalog #: 51809

Aliases

Interferon alpha-1/13; IFN-alpha-1/13; Interferon alpha-D; LeIF D

Background

Gene Name: IFNA1

NCBI Gene Entry: [3439/3447](#)

UniProt Entry: [P01562](#)

Application Information

Molecular Weight: Predicted, 21 kDa; observed, 23 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human IFN alpha 1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of IFN alpha 1 expression in HEK293T (A), A549 (B), MCF7 (C), mouse lung (D), mouse spleen (E), rat testis (F), rat spleen (G) whole cell lysates.