

Catalog #: 51818

Aliases

GPR172A; PAR1; RFT3; Solute carrier family 52, riboflavin transporter, member 2; Porcine endogenous retrovirus A receptor 1; PERV-A receptor 1; Protein GPR172A; Riboflavin transporter 3; hRFT3

Background

Gene Name: SLC52A2

NCBI Gene Entry: [79581](#)

UniProt Entry: [Q9HAB3](#)

Application Information

Molecular Weight: Predicted, 45 kDa; observed, 46 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human GPR172A

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of GPR172A expression in HEK293T (A), MCF7 (B) whole cell lysates.