

## Catalog #: 51824

### Aliases

SKD1; VPS42; Vacuolar protein sorting-associated protein 4B; Cell migration-inducing gene 1 protein; Suppressor of K(+) transport growth defect 1; Protein SKD1

### Background

Gene Name: VPS4B

NCBI Gene Entry: [9525](#)

UniProt Entry: [O75351](#)

### Application Information

Molecular Weight: Predicted, 49 kDa; observed, 49 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human

Applications Tested: Western blotting (WB)

### Immunogen

A synthesized peptide derived from human VPS4B

### Isotype

Rabbit IgG

### Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

### Storage

Store at -20 °C for one year.

### Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

**Note:** This product is for research use only.

### Validation Data

---


#### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

#### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)



Western blotting analysis of VPS4B expression in HEK293T (A), A549 (B), H1688 (C) whole cell lysates.