

Anti-Glutathione Peroxidase Rabbit Polyclonal Antibody



Catalog #: 51825

Aliases

Glutathione peroxidase 1; GPx-1; GSHPx-1; Cellular glutathione peroxidase

Background

Gene Name: GPX1

NCBI Gene Entry: [2876](#)

UniProt Entry: [P07203](#)

Application Information

Molecular Weight: Predicted, 22 kDa; observed, 22 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Glutathione Peroxidase

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of Glutathione Peroxidase expression in mouse liver (A), rat liver (B) whole cell lysates.