

Catalog #: 51848

Aliases

GLI; Zinc finger protein GLI1; Glioma-associated oncogene; Oncogene GLI

Background

Gene Name: GLI1

NCBI Gene Entry: [2735](#)

UniProt Entry: [P08151](#)

Application Information

Molecular Weight: Predicted, 117 kDa; observed, 140 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human GLI1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of GLI1 expression in Hela (A), mouse brain (B), A549 (C) whole cell lysates.