

Catalog #: 51895

Aliases

HOX1C; Homeobox protein Hox-A5; Homeobox protein Hox-1C

Background

Gene Name: HOXA5

NCBI Gene Entry: [3202](#)

UniProt Entry: [P20719](#)

Application Information

Molecular Weight: Predicted, 29 kDa; observed, 43,55 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human HOXA5

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of HOXA5 expression in COS7 (A), H1792 (B), Panc1 (C), MEF (D), PC12 (E) whole cell lysates.