

Catalog #: 51898

Aliases

FAM35A; Protein FAM35A; FAM35B; Protein FAM35B

Background

Gene Name: FAM35A/FAM35B

NCBI Gene Entry: [54537](#)

UniProt Entry: [Q86V20/A6NDK8](#)

Application Information

Molecular Weight: Predicted, 93 kDa; observed, 94 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human FAM35A/B

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of FAM35A/B expression in H1792 (A), H9C2 (B), AML12 (C) whole cell lysates.