

Catalog #: 51932

Aliases

CMKBR5; C-C chemokine receptor type 5; C-C CKR-5; CC-CKR-5; CCR-5; CCR5; CHEMR13; HIV-1 fusion coreceptor; CD195

Background

Gene Name: CCR5

NCBI Gene Entry: [1234](#)

UniProt Entry: [P51681](#)

Application Information

Molecular Weight: Predicted, 40 kDa; observed, 41 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, bovine, monkey, rabbit

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human CCR5

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of CCR5 expression in K562 (A), MCF7 (B) whole cell lysates.