

Anti-Androgen Receptor Rabbit Polyclonal Antibody



Catalog #: 51939

Aliases

DHTR; NR3C4; Androgen receptor; Dihydrotestosterone receptor; Nuclear receptor subfamily 3 group C member 4

Background

Gene Name: AR

NCBI Gene Entry: [367](#)

UniProt Entry: [P10275](#)

Application Information

Molecular Weight: Predicted, 99,110 kDa; observed, 110 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, rat, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Androgen Receptor

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of Androgen Receptor expression in Hela (A), MCF7 (B), rat testis (C) whole cell lysates.