

Catalog #: 51962

Aliases

UBB; Ubiquitin B; Polyubiquitin B; Polyubiquitin-B; FLJ25987; MGC8385; Epididymis Secretory Protein Li 50; HEL-S-50; UBC; Ubiquitin C; Polyubiquitin-C; HMG20; RPS27A; Ribosomal Protein S27a; Ubiquitin Carboxyl Extension Protein 80; UBCEP80; Uba80; S27A; ES31; Ubiquitin-Ribosomal Protein ES31 Fusion Protein; UBCEP1; Ubiquitin And Ribosomal Protein S27a; Ubiquitin-40S Ribosomal Protein S27a; Epididymis Luminal Protein 112; 40S Ribosomal Protein S27a; Ubiquitin-CEP80; Ubiquitin C; HEL112; CEP80; UBA80; UBC; UBA52; Ubiquitin A-52 Residue Ribosomal Protein Fusion Product 1; CEP52; HUBCEP52; RPL40; L40; Ubiquitin-Ribosomal Protein EL40 Fusion Protein; Ubiquitin Carboxyl Extension Protein 52; Ubiquitin-52 Amino Acid Fusion Protein; Ubiquitin-CEP52; MGC126879; MGC126881; MGC57125; Ubiquitin-60S Ribosomal Protein L40; 60S Ribosomal Protein L40; Ribosomal Protein L40; UBCEP2

Background

Gene Name: UBA52/RPS27A/UBB/UBC

NCBI Gene Entry: [7311/7314/7316/6233](#)

UniProt Entry: [P0CG47/P0CG48/P62979/P62987](#)

Application Information

Molecular Weight: Predicted, 25,77,17,14 kDa; observed, 43 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Ubiquitin

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-Ubiquitin antibody (Cat#51962). Total lysates (30 µg) were loaded and separated by SDS-PAGE. The blot was incubated with anti-Ubiquitin antibody (Cat#51962, 1:2,500) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM