

Catalog #: 51967

Aliases

VEGF; Vascular endothelial growth factor A; VEGF-A; Vascular permeability factor; VPF

Background

Gene Name: VEGFA

NCBI Gene Entry: [7422](#)

UniProt Entry: [P15692-4](#)

Application Information

Molecular Weight: Predicted, 27 kDa; observed, 21,42 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, pig

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human VEGF165

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of VEGF165 expression in mouse embryo (A), rat brain (B), LO2 (C), Beas2B (D), A549 (E), Hela (F), SGC7901 (G) whole cell lysates.