

Catalog #: 51975

Aliases

CABP9K; CALB3; S100D; Protein S100-G; Calbindin-D9k; S100 calcium-binding protein G; Vitamin D-dependent calcium-binding protein intestinal; CABP

Background

Gene Name: S100G

NCBI Gene Entry: [795](#)

UniProt Entry: [P29377](#)

Application Information

Molecular Weight: Predicted, 9 kDa; observed, 16 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, pig

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human S100-G

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of S100-G expression in rat spleen (A) whole cell lysates.