

Catalog #: 52056

Aliases

PRR3; Poliovirus receptor-related protein 3; CDw113; Nectin-3; CD113

Background

Gene Name: PVRL3

NCBI Gene Entry: [25945](#)

UniProt Entry: [Q9NQS3](#)

Application Information

Molecular Weight: Predicted, 61 kDa; observed, 95 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human CD113

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of CD113 expression in Hela (A), A549 (B), C6 (C), AML12 (D) whole cell lysates.