

Catalog #: 52087

Aliases

Interleukin-7 receptor subunit alpha; IL-7 receptor subunit alpha; IL-7R subunit alpha; IL-7R-alpha; IL-7RA; CDw127; CD127

Background

Gene Name: IL7R

NCBI Gene Entry: [3575](#)

UniProt Entry: [P16871](#)

Application Information

Molecular Weight: Predicted, 51 kDa; observed, 68,90 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human CD127

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of CD127 expression in SKOV3 (A), HEK293T (B), rat kidney (C), mouse testis (D) whole cell lysates.



Western blotting analysis of CD127 expression in A549 (A), SGC7901 (B), C6 (C), CT26 (D) whole cell lysates.