

# Anti-STAT3 (Acetyl-K87) Rabbit Polyclonal Antibody



**Catalog #: 52089**

## Aliases

APRF; Signal transducer and activator of transcription 3; Acute-phase response factor

## Background

Gene Name: STAT3

NCBI Gene Entry: [6774](#)

UniProt Entry: [P40763](#)

## Application Information

Molecular Weight: Predicted, 88 kDa; observed, 88 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, chicken, pig

Applications Tested: Western blotting (WB)

## Immunogen

A synthesized peptide derived from human STAT3 protein

## Isotype

Rabbit IgG

## Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

## Storage

Store at -20 °C for one year.

## Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

**Note:** This product is for research use only.

## Validation Data

---


### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)



Western blotting analysis of STAT3 (Acetyl-K87) expression in Hela (A), H446 (B), mouse lung (C), mouse kidney (D), rat lung (E), rat kidney (F) whole cell lysates.