

Catalog #: 52111

Aliases

c Myc; c Myc Epitope Tag; c Myc Tag; Myc Epitope Tag; Myc proto-oncogene protein; Transcription factor p64

Background

Gene Name: N/A

UniProt Entry: [N/A](#)

Application Information

Molecular Weight: B.

Clonality: Rabbit polyclonal antibody

Species Reactivity: N/A

Applications Tested: Western blotting (WB), immunoprecipitation (IP)

Immunogen

A synthesized peptide derived from EQKLISEEDL.

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:2,000

Immunoprecipitation (IP): 1:10-1:100

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of Anti-Myc-tag Antibody against HEK293T cells transfected with vector overexpressing Myc tag (A) and untransfected (B).