

Catalog #: 52115

Aliases

Mitogen-activated protein kinase 6; PROSITE-ProRule annotation; AtMPK6; MAP kinase 6

Background

Gene Name: MPK6

NCBI Gene Entry: [818982](#)

UniProt Entry: [Q39026](#)

Application Information

Molecular Weight: Predicted, 45 kDa; observed, 45 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: At

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from the N-term region of arabidopsis thaliana MPK6

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of MPK6 expression in Arabidopsis thaliana (A) whole cell lysates.