

Catalog #: 52116

Aliases

Serine/threonine-protein phosphatase PP1-beta catalytic subunit; PP-1B; PPP1CD

Background

Gene Name: PPP1CB

NCBI Gene Entry: [5500](#)

UniProt Entry: [P62140](#)

Application Information

Molecular Weight: Predicted, 37 kDa; observed, 37 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human PP1 beta

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of PP1 beta expression in HL60 (A), PC3 (B), Hela (C), mouse heart (D) whole cell lysates.