

Catalog #: 52120

Aliases

PRSS10; Transmembrane protease serine 2; Serine protease 10

Background

Gene Name: TMPRSS2

NCBI Gene Entry: [7113](#)

UniProt Entry: [O15393](#)

Application Information

Molecular Weight: Predicted, 53 kDa; observed, 65,54,32 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human TMPRSS2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of TMPRSS2 expression in HEK293T (A), K562 (B), MCF7 (C), mouse spleen (D), rat kidney (E), rat spleen (F) whole cell lysates.