

Catalog #: 52129

Aliases

ATGL; Patatin-like phospholipase domain-containing protein 2; Adipose triglyceride lipase; Calcium-independent phospholipase A2; Desnutrin; IPLA2-zeta; Pigment epithelium-derived factor; TTS2.2; Transport-secretion protein 2; TTS2

Background

Gene Name: PNPLA2

NCBI Gene Entry: [57104](#)

UniProt Entry: [Q96AD5](#)

Application Information

Molecular Weight: Predicted, 55 kDa; observed, 55,47 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human ATGL

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of ATGL expression in A549 (A), HepG2 (B) whole cell lysates.