

Catalog #: 52144

Aliases

G11; RP1; Serine/threonine-protein kinase 19; Protein G11; Protein RP1

Background

Gene Name: STK19

NCBI Gene Entry: [8859](#)

UniProt Entry: [P49842](#)

Application Information

Molecular Weight: Predicted, 40 kDa; observed, 40,28 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human STK19

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of STK19 expression in Jurkat (A), HepG2 (B), mouse liver (C), rat lung (D), rat brain (E) whole cell lysates.