

Catalog #: 52155

Aliases

CLEC13C; Secretory phospholipase A2 receptor; PLA2-R; PLA2R; 180 kDa secretory phospholipase A2 receptor; C-type lectin domain family 13 member C; M-type receptor

Background

Gene Name: PLA2R1

NCBI Gene Entry: [22925](#)

UniProt Entry: [Q13018](#)

Application Information

Molecular Weight: Predicted, 152,168 kDa; observed, 168 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human PLA2R

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of PLA2R expression in SHSY5Y (A), mouse liver (B), mouse heart (C) whole cell lysates.