

Catalog #: 52177

Aliases

CCDC137; Coiled-Coil Domain Containing 137; Coiled-Coil Domain-Containing Protein 137; MGC16597; Hepatocellular Carcinoma Related Protein; Retinoic Acid Resistance Factor; CPERP-B; RaRF

Background

Gene Name: CCDC137

NCBI Gene Entry: [339230](#)

UniProt Entry: [Q6PK04](#)

Application Information

Molecular Weight: Predicted, 33 kDa; observed, 32 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human CC137

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

