

Catalog #: 52281

Aliases

ALS4; AOA2; bA479K20.2; SCAR1

Background

Gene Name: SETX

NCBI Gene Entry: [23064](#)

UniProt Entry: [Q7Z333](#)

Application Information

Molecular Weight: Predicted, 35 kDa; observed, 48 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human

Applications Tested: Western blotting (WB), flow cytometry (FCM), immunocytochemistry (IC), immunoprecipitation(IP), immunohistochemistry (IHC)

Immunogen

A synthesized peptide derived from human Senataxin protein

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Flow Cytometry (FCM): 1:10-1:1,000

Immunocytochemistry (IC): 1:100

Immunoprecipitation (IP): 1:500,

Western Blotting (WB): 1:500

Immunohistochemistry (IHC): 1:200-1:500

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Validation Data

