

Catalog #: 52283

Aliases

GPR17; G Protein-Coupled Receptor 17; Uracil Nucleotide/Cysteinyl Leukotriene Receptor; UDP/CysLT Receptor; P2Y-Like Receptor; R12; G-Protein Coupled Receptor 17

Background

Gene Name: GPR17

NCBI Gene Entry: [2840](#)

UniProt Entry: [Q13304](#)

Application Information

Molecular Weight: Predicted, 41 kDa

Clonality: Rabbit polyclonal antibody

Clone ID: 24GB8620

Species Reactivity: Human

Applications Tested: Immunohistochemistry (IHC)

Immunogen

A synthesized peptide derived from human GPR17

Isotype

Mouse IgG2b

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Immunohistochemistry(IHC): 1:150~600

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

