

Catalog #: 52287

Aliases

RBPJL; Recombination Signal Binding Protein For Immunoglobulin Kappa J Region Like; RBPSUHL; SUHL; Recombining Binding Protein Suppressor Of Hairless-Like Protein; Transcription Factor RBP-L; RBP-L; RBPL; SUH; Recombination Signal Binding Protein For Immunoglobulin Kappa J Region-Like; Recombining Binding Protein Suppressor Of Hairless (Drosophila)-Like

Background

Gene Name: RBPJL

NCBI Gene Entry: [11317](#)

UniProt Entry: [Q9UBG7](#)

Application Information

Molecular Weight: Predicted, 57 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human RBPJL

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

