

Anti-IL1RL2 (spesolimab biosimilar) Monoclonal Antibody



Catalog #: 52295

Aliases

IL1RL2; Interleukin 1 Receptor Like 2; IL1R-Rp2; IL1RRP2; IL-36R; Interleukin-1 Receptor-Related Protein 2; Interleukin-1 Receptor-Like 2; IL-36 Receptor; IL-1Rrp2; IL-1 Receptor-Related Protein 2; EC 3.2.2.6

Background

Gene Name: IL1RL2

NCBI Gene Entry: 8808

UniProt Entry: Q9HB29

Application Information

Molecular Weight: Predicted, 65 kDa

Clonality: monoclonal antibody

Species Reactivity: Human

Applications Tested: Flow cytometry (FCM)

Immunogen

A synthesized peptide derived from human IL1RL2

Isotype

Human IgG1

Storage Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5% – 8% trehalose is added as protectants before lyophilization.

Storage

Store at -20°C to -80°C for 12 months in lyophilized form.

Recommended Dilutions

Flow Cytometry (FCM): 1:100

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Flow cytometry analysis with 1 μ g/mL Anti-IL1RL2 (spesolimab biosimilar) mAb on Expi293 cells transfected with Human IL1RL2 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).