

Invivo Anti-Mouse Ly6G Rat Monoclonal Antibody



Catalog #: 52302

Aliases

Lymphocyte antigen 6G; Ly-6G; Ly-6G.1

Background

Gene Name: Ly6G

NCBI Gene Entry: [10090](#)

UniProt Entry: [P35461](#)

Application Information

Molecular Weight: Predicted, 14 kDa

Clonality: Rat monoclonal antibody

Clone ID: 24GB13900

Species Reactivity: Mouse

Applications Tested: In vivo neutrophil depletion, in vivo MDSC depletion

Isotype

Rat IgG2a,k

Storage Buffer

Supplied in PBS (pH 7.4) containing no preservativ.

Storage

12 months from date of receipt / reconstitution, 2 to 8 °C as supplied

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

