

Catalog #: 52303

Aliases

RBM4, RNA-binding protein 4, Lark homolog, hLark, RNA-binding motif protein 4, RNA-binding motif protein 4a, RBM4, RBM4A

Background

Gene Name: RBM4

NCBI Gene Entry: [5936](#)

UniProt Entry: [Q9BWF3](#)

Application Information

Molecular Weight: Predicted, 40 kDa; observed, 37 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human RBM4

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-RBM4 antibody (Cat#52303). Total lysates (30 μ g) were loaded and separated by SDS-PAGE. The blot was incubated with anti-RBM4 antibody (Cat#52303, 1:2,500) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQTM ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM