

Catalog #: 52323

Aliases

IDO2; Indoleamine 2,3-Dioxygenase 2; INDOL1; Indoleamine-Pyrrole 2,3-Dioxygenase-Like Protein 1; Indoleamine-Pyrrole 2,3 Dioxygenase-Like 1; Indoleamine 2,3-Dioxygenase-Like Protein 1; Indoleamine 2,3-Dioxygenase-Like 1 Protein; EC 1.13.11.-; IDO-2

Background

Gene Name: IDO2

NCBI Gene Entry: [169355](#)

UniProt Entry: [Q6ZQW0](#)

Application Information

Molecular Weight: Predicted, 45 kDa; observed, 45 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human IDO2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of HeLa, HepG2 cells using anti-INDOL1 antibody (Cat#52323). The blot was incubated with anti-INDOL1 antibody (Cat#52323) and HRP-conjugated goat anti-rabbit secondary antibody (1:20,000) respectively.