

Catalog #: 52353

Aliases

LCAT; Lecithin-Cholesterol Acyltransferase; 1-Alkyl-2-Acetylgllycerophosphocholine Esterase; Phosphatidylcholine--Sterol O-Acyltransferase; Phosphatidylcholine-Sterol Acyltransferase; Platelet-Activating Factor Acetylhydrolase; Phospholipid-Cholesterol Acyltransferase; PAF Acetylhydrolase; EC 2.3.1.43; Testicular Secretory Protein Li 24; EC 3.1.1.47

Background

Gene Name: LCAT

NCBI Gene Entry: [3931](#)

UniProt Entry: [P04180](#)

Application Information

Molecular Weight: Predicted, 50 kDa; observed, 55 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB16055

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human LCAT

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

