

Catalog #: 52354

Aliases

PPARG; Peroxisome Proliferator Activated Receptor Gamma; NR1C3; Peroxisome Proliferator-Activated Receptor Gamma; PPARgamma; PPARG1; PPARG2; Nuclear Receptor Subfamily 1 Group C Member 3; PPAR-Gamma; Peroxisome Proliferator-Activated Receptor-Gamma Splicing; Peroxisome Proliferative Activated Receptor, Gamma; Peroxisome Proliferator-Activated Receptor-Gamma 5; PPARG5; CIMT1; FPLD3; GLM1

Background

Gene Name: PPARG

NCBI Gene Entry: [5468](#)

UniProt Entry: [P37231](#)

Application Information

Molecular Weight: Predicted, 53-57 kDa; observed, 53-57 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 24GB16060

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), immunohistochemistry (IHC)

Immunogen

A synthesized peptide derived from human PPAR γ

Isotype

Mouse IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:2,000

Immunohistochemistry (IHC): 1:50-1:200

Note: This product is for research use only.

Validation Data

