

Catalog #: 52357

Aliases

IL17RA; Interleukin 17 Receptor A; IL-17RA; CDw217; HIL-17R; CD217; IL17R;
Interleukin-17 Receptor A; IL-17 Receptor A; Interleukin 17 Receptor; CD217 Antigen;
CANDF5; IMD51

Background

NCBI Gene Entry: [23765](#)

UniProt Entry: [Q96F46](#)

Application Information

Molecular Weight: Predicted, 96 kDa, observed, 90 kDa

Clonality: Rabbit polyclonal antibody

Clone ID: N/A

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human IL17RA

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,500

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-IL17RA antibody (Cat#52357). Total lysates (30 µg) were loaded and separated by SDS-PAGE. The blot was incubated with anti-IL17RA antibody (Cat#52357, 1:2,500) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX:+1-540-855-7041

WWW.GENUINBIOTECH.COM