

Catalog #: 52372

Aliases

FLNB; Filamin B; ABP-278; TABP; TAP; FLN1L; FH1; Larsen Syndrome 1 (Autosomal Dominant); Actin-Binding-Like Protein; Actin Binding Protein 278; Thyroid Autoantigen; Filamin Homolog 1; Filamin B, Beta; ABP-280 Homolog; Beta-Filamin; Filamin-B; Filamin-3; FLN-B; LRS1; Filamin B, Beta (Actin Binding Protein 278); Truncated Actin-Binding Protein; Truncated ABP; Beta Filamin; ABP-280; FLN3; AOI; SCT; Fh1

Background

Gene Name: FLNB

NCBI Gene Entry: [2317](#)

UniProt Entry: [O75369](#)

Application Information

Molecular Weight: Predicted, 278 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human FLNB

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

