

Catalog #: 52373

Aliases

MYH9; Myosin Heavy Chain 9; NMHC-II-A; NMMHCA; EPSTS; FTNS; MHA; Cellular Myosin Heavy Chain, Type A; Nonmuscle Myosin Heavy Chain II-A; Non-Muscle Myosin Heavy Chain IIA; Non-Muscle Myosin Heavy Chain A; NMMHC-IIA; Myosin-9; DFNA17; Myosin, Heavy Polypeptide 9, Non-Muscle; Non-Muscle Myosin Heavy Polypeptide 9; Myosin Heavy Chain, Non-Muscle IIA; Myosin, Heavy Chain 9, Non-Muscle; Non-Muscle Myosin Heavy Chain 9; Nonmuscle Myosin IIA2; NMMHC II-A; NMMHC-A; BDPLT6; MATINS

Background

Gene Name: MYH9

NCBI Gene Entry: [4627](#)

UniProt Entry: [P35579](#)

Application Information

Molecular Weight: Predicted, 227 kDa; observed, 215 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human MYH9

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Validation Data

